

enbio


# Global presence



Enbio Group AG is a part  
of global investment group  
with a turnover of EUR 200M.



B, like the best in class



**Vacuum  
(Class B)**

Class S

Class N

# Autoclave classes, explained

	B	S	N
Drying	✓	-	-
Pouches	✓	-	-
<u>Any</u> instrument type	✓	instrument mfr. approved only	-
Pre-vacuums	3	1	-



Dentistry



Veterinary




Podiatry



Surgery



Piercing and Tattoo



Beauty

# The fastest

Why would you waste an hour on something you can do in 7 minutes? Enbio B Class air-removal (pre-vacuum) autoclave is the fastest sterilizer on the market. It is also gloriously intuitive and easy to use.

Sterilization process takes **4 min.**

Sterilization between visits.


Immediately ready for the next sterilization.



# Spacious

It will fit any interior thanks to its compact size (56x25x16 cm). At the same time, it has the space for all the standard instruments, in their packaging or without it.

## Enbio




Load space  
28x18x3 cm



x 8

## Enbio Pro



Load space  
28x18x7 cm



x 16







	Typical autoclave	Enbio PRO	Enbio
Tool pouches / h	50	50	32
Vacuum / drying	✓	✓	✓

# Silent

It generates a low noise level of 40 dB, which allows for comfortable work nearby the device.

Sterilization while being with your patient.



Mobile



enbio.

Designed   
in Switzerland



reddot design award  
winner



enbio.

# Enbio

Fastest autoclave in the world for instruments sterilization  
in medical and beauty sector.

Unwrapped FAST / 134°C	Wrapped Class B / 134°C	Wrapped Class B / 121°C
<b>7 min</b>	<b>15 min</b>	<b>30 min</b>



# Enbio Pro

Fastest autoclave in the world for instruments sterilization  
in medical and beauty sector.

Unwrapped  
FAST / 134°C

**9 min**

Wrapped  
Class B / 134°C

**18 min**

Wrapped  
Class B / 121°C

**30 min**

Wrapped / PRION  
Class B / 134°C

**43 min**





One touch – and ready!





# Features

## Power supply



Enbio works with any typical electrical outlet.

## Water supply



Enbio uses distilled water from external tank. With Enbio Magic Filter, water can be recycled.

## Data archiving



All process data are automatically written to a USB disk.

## Process reports



Enbio Data Viewer (included) allows generating and printing detailed process reports from the data on USB.


# 80 ml water per process

(Vacuum / class B)



enbio.

# Enbio Magic Filter: forget about water



**1 liter  
of demineralized  
water**

**1000 processes**




Designed  
in Switzerland



FDA  
Cleared



# Enbio Care



Perfect technology meets perfect care. All Enbio devices are covered with world's unique warranty and care program. We offer a special qr troubleshooting system for our users.

# Marketing support

- Regular product trainings
- Marketing kit:
  - photostock
  - brochures
  - leaflets
- Marketing support with expert



enbio.

# Field-proven product

36,000 happy users in just 6 years



enbio.



Thank you